


Dentapreg™ Implant Overdenture Support

English

- 1) Insert the flow composite to the post on area approximately 3x3 mm.




- 2) Remove the paper from the strip and the end of the strip insert to the layer of the flow composite on the post - hold it with pincers through the foil and cure it pointwise only for 3-5 seconds.



- 3) Wind around the another post with fiber strip.




- 4) Press the both parts of the strip with pincers together through the foil and pointwise light cure it (from 3 to 5 seconds).



- 5) Finish the wind fiber strip around the second and the first post so that the end of the strip is ended at the half of the distance between the posts.



- 6) Press the parts of the strip with the pincers together through the foil and light cure it (1 cm for 40 seconds).



- 7) Cover the end of the strip with a thin layer of the flow composite. Cure the flow composite according to manufacturer instructions.



- 8) Apply a flow composite on fiber composite closes to the post. Remove the paper from SFM strip. The beginning of the SFM strip insert to the flow composite and pointwise light cure it.



- 9) Wind SFM strip around the FRC construction.



- 10) Press the winded SFM strip to the FRC construction and light cure it.



- 11) Cover the end of the SFM strip with a thin layer of flow and light cure it according to manufacturer instructions. The steps 8, 9, 10 repeat as necessary (length of construction) however wind the SFM strip from the second post to the center of FRC construction.




- 12) Fill the cavity between the strip and posts with flow composite and light cure it according to manufacturer instructions.



Dentapreg™ Implant Overdenture Support

English

- 1) Insert the flow composite to the post on area approximately 3x3 mm.




- 2) Remove the paper from the strip and the end of the strip insert to the layer of the flow composite on the post - hold it with pincers through the foil and cure it pointwise only for 3-5 seconds.



- 3) Wind around the another post with fiber strip.



- 4) Press the both parts of the strip with pincers together through the foil and pointwise light cure it (from 3 to 5 seconds).



- 5) Finish the wind fiber strip around the second and the first post so that the end of the strip is ended at the half of the distance between the posts.



- 6) Press the parts of the strip with the pincers together through the foil and light cure it (1 cm for 40 seconds).



- 7) Cover the end of the strip with a thin layer of the flow composite. Cure the flow composite according to manufacturer instructions.



- 8) Apply a flow composite on fiber composite closes to the post. Remove the paper from SFM strip. The beginning of the SFM strip insert to the flow composite and pointwise light cure it.



- 9) Wind SFM strip around the FRC construction.



- 10) Press the winded SFM strip to the FRC construction and light cure it.



- 11) Cover the end of the SFM strip with a thin layer of flow and light cure it according to manufacturer instructions. The steps 8, 9, 10 repeat as necessary (length of construction) however wind the SFM strip from the second post to the center of FRC construction.



- 12) Fill the cavity between the strip and posts with flow composite and light cure it according to manufacturer instructions.



Dentapreg™ Vyztužení náhrady na implantatech

Dentapreg™ Vyztužení náhrady na implantatech

Česky



1) Aplikujte flow kompozit na první čep.



(Na průčelu přibízne 3 x 3 mm).



3) Dlouhého časného druhého čepu.



4) Sískněte čast páska k sobě (pres foil) a zároveň ho bodově vytrhněte.



5) Dlouhého časného páska k sobě a zároveň vytrhněte.



6) Sískněte konc páska ke zbytku páska a zároveň vytrhněte.



7) Konc páska překývají slabou vstavu flow kompozitu a lehce vytrhněte v souladu s návodem vyrobce.



8) Aplikujte flow kompozit na páska co nejjízde a bodově ho vytrhněte.



9) Odložte SFM páska kolem FRC konstrukce směrem k jeřimu středu.



10) Po odložení páska ho přilákejte přezutou k FRC konstrukci a vytrhněte ho.



11) Konc obloženého SFM páska překývají slabou vstavu flow kompozitu a vytvořte ho v souladu s doporučením vyrobce.



12) Zapište kavity mezi páska a čepy pomocí flow kompozitu a vytvořte ho v souladu s doporučením vyrobce.



13) Konc obloženého SFM páska překývají slabou vstavu flow kompozitu a vytvořte ho v souladu s doporučením vyrobce.



14) Konc obloženého SFM páska překývají slabou vstavu flow kompozitu a vytvořte ho v souladu s doporučením vyrobce.



15) Aplikujte flow kompozit na páska co nejjízde a bodově ho vytrhněte.



16) Sískněte konc páska ke zbytku páska a zároveň vytrhněte.



17) Konc páska překývají slabou vstavu flow kompozitu a lehce vytrhněte v souladu s návodem vyrobce.



18) Aplikujte flow kompozit na páska co nejjízde a bodově ho vytrhněte.



19) Odložte obloženého páska kolem druhého čepu a zároveň vytrhněte.



20) Po odložení páska ho přilákejte přezutou k FRC konstrukci a vytrhněte ho.



21) Zapište kavity mezi páska a čepy pomocí flow kompozitu a vytvořte ho v souladu s doporučením vyrobce.



22) Sískněte silikonový papír z jednoho konc vytvořte po dobu 3-5 s.



23) Odložte páska kolem druhého čepu.



24) Aplikujte flow kompozit na první čep.



25) Dokončete obloženého páska kolem druhého čepu a zároveň vytrhněte ho tak aby jeho koncové zasobování do středu vzdleností mezi a prvního čepu a síťnicí nebo tak aby jeho koncové zasobování do středu vzdleností mezi oběma čepy.

26) Sískněte konc pásku ke zbytku páska a zároveň vytrhněte.

27) Konc páska překývají slabou vstavu flow kompozitu a lehce vytrhněte v souladu s návodem vyrobce.

28) Aplikujte flow kompozit na páska co nejjízde a bodově ho vytrhněte.

29) Odložte obloženého páska kolem druhého čepu a zároveň vytrhněte.

30) Sískněte konc pásku ke zbytku páska a zároveň vytrhněte.

31) Dokončete obloženého páska kolem druhého čepu a zároveň vytrhněte ho tak aby jeho koncové zasobování do středu vzdleností mezi a prvního čepu a síťnicí nebo tak aby jeho koncové zasobování do středu vzdleností mezi oběma čepy.

32) Sískněte silikonový papír z jednoho konc vytvořte po dobu 3-5 s.

33) Odložte páska kolem druhého čepu.

34) Aplikujte flow kompozit na první čep.

35) Na průčelu přibízne 3 x 3 mm.

36) Česky

navod implantaty.cz.indd 1
nařízení 15.10.2009 14:15:45